


# STERILE SPECIMEN COLLECTION SWAB

Nasopharyngeal Swab with 80mm Breakpoint



As we navigate through these unprecedented times, Applied Medical is committed to rapidly responding to the needs of the healthcare community. We are all in this together, fighting against COVID-19.

# NASOPHARYNGEAL SWAB WITH 80mm BREAKPOINT

Applied Medical's sterile specimen collection swab is indicated for the collection of nasopharyngeal clinical samples to support COVID-19 testing. The swab consists of a medical-grade polystyrene handle and polyester tip.

	Applied Medical Nasopharyngeal Swab (NPS21)	Traditional Flocked Swab
Tip Material	Polyester	Polyester
Tip Diameter	3.5mm (Compressed: 0.85mm)	3.56mm
Breakpoint	80mm	80mm
Shaft Material	Polystyrene	Polystyrene
Length	6.0"	6.0"


## FEATURES

- Sterile, individually packaged swabs
- Made in the USA with materials sourced domestically

## ORDER INFORMATION

Model #	Description	Qty/Box
NPS21	Sterile specimen collection swab	200**

\*\*Controlled by weight



Applied Medical is pleased to offer nasopharyngeal swabs direct to U.S. healthcare facilities and government agencies. For more information or to inquire about placing an order, please contact [COVIDResponse@appliedmedical.com](mailto:COVIDResponse@appliedmedical.com).

\*Compared to Traditional Flocked swabs

© 2020 Applied Medical Resources Corporation. All rights reserved. Applied Medical, the Applied Medical logo design and marks designated with a ® are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States, and/or the European Union. 327687-EN-USA-A

